


Solutions for Evaluation of CAR Expression

Specific detection of CAR expression using target antigens



Content

P 01	I. Introduction
P 01	II. Evaluation of CAR Expression
P 01	III. CAR Detection Strategy and Product Design
P 02	• Three Detection Methods
P 02	• FDA DMF Filed Proteins
P 03	IV. Case Studies
P 07	V. Product List of CAR-T Targets
P 07	• Fluorescent-labeled Proteins
	◆ FITC-labeled proteins
	◆ PE-labeled proteins
P 08	• Biotinylated Proteins
P 10	• Unconjugated Proteins

CD22

ROR1

MSLN

GPC3

BCMA


EGFRVIII

CD19

I. Introduction

The chimeric antigen receptor T (CAR-T) cell therapy is a new treatment for a variety of cancers. The idea is to take out the T-cells from the patient, and genetically engineer the cells to make them express a chimeric antigen receptor (CAR) recognizing a specific tumor-associated antigen (TAAs). As a result, the CAR-expressing T cells, when reintroduced into the patient's body, will target and eliminate the TAA-expressing tumor cells.

ACROBiosystems has developed an extensive collection of recombinant proteins to support CAR-T therapy development. This growing list of proteins includes many fluorescent-labeled target antigens and pre-biotinylated proteins that are uniquely suitable for evaluation of CAR expression. In addition, we also supply hard-to-make proteins such as BCMA, CD19, ROR1, and EGFRVIII.



II. Evaluation of CAR Expression

Evaluating CAR expression is an essential step in the production of CAR-T cells. This is often done by flow cytometry, using protein L, anti-Fab antibodies, anti-idiotypic Antibodies or target antigens as detection antibodies. Among these common choices, target antigens are widely considered to be the best option, because it offers high specificity and minimal background staining.

Reagents	Mechanism	Pros	Cons
Target Antigens, Anti-idiotypic Antibodies	Specifically bind to the antigen-binding domains of CARs.	High specificity; Minimal background staining.	Each unique CAR has to be stained with corresponding antigens.
Protein L	Binds to the kappa light chain of immunoglobulin.	Universal.	High background staining; Cannot detect the anti-lambda light chain CAR.
Anti-Fab antibody	Binds to the Fab portion of immunoglobulin.	Universal.	High background staining.

III. CAR Detection Strategy and Product Design

Currently, target antigens for CAR-T cells are most widely used to determine the expression of CARs on gene-modified lymphocytes by flow cytometry. The limitations of these reagents are that many of them are not commercially available. In an effort to fulfill these needs, we have developed an extensive collection of CAR-T target antigens includes many fluorescent-labeled proteins and pre-biotinylated proteins that are uniquely suitable for evaluation of CAR expression by flow cytometry.


Detection Methods

■ Direct detection

- Target antigens are pre-labeled with fluorescent dye.
- Processing time can be reduced by the use of direct-labeled proteins.
- Non-specific reaction of a secondary antibody is eliminated.

ACRO's specially designed products:

FITC-labeled proteins;
PE-labeled proteins;




■ Biotin-streptavidin based detection

- Target antigens are pre-labeled with biotin and detected by labeled streptavidin (the biotin-avidin complex).
- Streptavidin labeled with fluorochromes can bind biotinylated proteins with a high degree of affinity and specificity, amplifying the signal and improving the detection sensitivity and specificity.

ACRO's specially designed products:

Avitag™ biotinylated proteins;
Chemically biotinylated proteins;




■ Indirect detection

- Target antigens are designed to carry a specific tag and detected using a secondary antibody (anti-epitope tag antibody) labeled with a fluorophore.
- Non-specific reaction of a secondary antibody may occur.

ACRO's specially designed products:

Fc-tag fusion proteins;
His-tag fusion proteins;



FDA DMF Filed Proteins

ACROBiosystems has submitted DMF for its recombinant CD19 and BCMA proteins to FDA, and filed the DMF number as 034936. You can cite this DMF number to shorten the time of preparing the filing of documents to support your IND or BLA.

■ How to get our DMF authorization

If our DMF filed proteins have been used in your drug development process, you can request that we provide DMF authorization to FDA in support of a submission or filing that you have made to the FDA. Please submit your request to ACROBiosystems by leaving your information here. [Request for Authorization](#)

* If you have FDA DMF filing request for other products, please [contact us](#).

IV. Case Studies

Evaluation of Anti-CD19 CAR Expression with FITC-labeled CD19

■ Reagents

FITC-labeled Human CD19 (20-291) Protein (ACROBiosystems, Cat. No. [CD9-HF2H2](#));

■ Samples


R1013-C6 cells (Transfected 293 cells expressing the anti-CD19 [FCM63] scFv & RFP tag).

■ Protocol

1. Culture R1013-C6 cells in DMEM medium with 10% FBS in the CO₂ incubator (at 37°C, 5% CO₂).
2. Harvest the cells and wash the cells once by wash buffer.
3. Count the cells number and the viability, aliquot up 2e5 live cells (Anti-CD19-scFv positive cell is 98%) into each tube. (Note: the cell viability must be ≥ 95%.)
4. Add 100 μl, 10 μg/ml of FITC-labeled Human CD19 (20-291) Protein or FITC-labeled Protein control into each tube, incubating at 4°C for 1 hour.
5. Wash the cells 3 times by wash buffer and resuspend the cells in 200 μl PBS per sample.
6. Transfer the cells into flow tube and detect by Flow cytometry.
7. Analyze result using FACS Celesta software and FCS Express 6 Flow software.

■ Results

The data showed that the expression level of anti-CD19 scFv on the surface of R1013-C6 cells was 95.37%.



293 cells were transfected with anti-CD19-scFv and RFP tag. 2e5 of the cells were stained with B. FITC-labeled Human CD19 (20-291) (Cat. No. [CD9-HF2H2](#), 10 μg/ml) and C. FITC-labeled protein control. A. Non-transfected 293 cells and C. FITC-labeled protein control were used as negative control. RFP was used to evaluate CAR (anti-CD19-scFv) expression and FITC was used to evaluate the binding activity of FITC-labeled Human CD19 (20-291) (Cat. No. [CD9-HF2H2](#)).

The protocol can be provided and feel free to send email to cart@acrobiosystems.com to request.

[Click here](#)

Evaluation of Anti-BCMA CAR Expression with Biotinylated BCMA

■ Reagents

Biotinylated human BCMA protein, Fc & Avi Tag (ACROBiosystems, Cat. No. [BC7-H82F0](#));

■ Samples


Anti-BCMA CAR-transduced human primary T-cells.

■ Brief Protocol

1. Human T cells transfected with anti-BCMA CAR were harvested 3 days after the transfection.
2. Aliquot up to 1e6 cells into centrifuge tube and wash the cells twice with FACS buffer.
3. Resuspend cells in 50 μ l of diluted biotinylated human BCMA (ACROBiosystems, Cat. No. BC7-H82F0) (prepared in FACS buffer at 8 μ g/ml) and incubate at 4°C for 30 minutes.
4. Wash cells twice with FACS buffer.
5. Resuspend cells in 50 μ l of diluted PE Streptavidin (Biolegend, Cat. No. 405204) (prepared in FACS buffer at 1:50 dilution) and incubate at 4°C for 30 minutes in the dark.
6. Wash cells twice with FACS buffer and resuspend the cells in 400 μ l PBS.
7. Transfer the cells into flow tube and analyze on BD FACSCalibur™ flow cytometer using FCS Express 6 Plus software.

■ Results

The data showed that the expression of anti-BCMA CARs on transduced T cell surface from donor 1 and donor 2 were 52.72% and 73.49%, respectively.



Human T cells were transfected with anti-BCMA CAR and cultured for 3 days. Three days post-transfection, 1e6 cells were first incubated with 50 μ l biotinylated human BCMA protein (Cat. No. [BC7-H82F0](#), 8 μ g/ml), washed and then stained with PE Streptavidin. (Data are kindly provided by PREGENE Biopharma)

The protocol can be provided and feel free to send email to cart@acrobiosystems.com to request.

[Click here](#)

Evaluation of Anti-CD19 CAR Expression with FITC-labeled Anti-FMC63 scFv Antibody

■ Reagents

FITC-Labeled Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (ACROBiosystems, Cat. No. [FM3-FY45](#)).

■ Samples


Anti-CD19 CAR-293 cells

■ Brief Protocol

1. Culture Anti-CD19 CAR-293 cells in DMEM medium with 10% FBS in the CO₂ incubator (at 37 °C, 5% CO₂).
2. Harvest the cells and wash the cells once by FACS buffer.
3. Count the cells number and the viability, aliquot up 2×10^5 live cells into each tube. (Note: the cell viability must $\geq 95\%$.)
4. Dilute FITC-Labeled Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (ACROBiosystems, Cat. No. FM3-FY45) in FACS buffer to get the working solution just before the assay, and then add 100 μ L of the working solution into the tube with cell pellet. Mix well and incubate at 4°C for 60 minutes.
5. Wash the cells 3 times by FACS buffer and resuspend the cell pellet in 200 μ L PBS per sample.
6. Transfer the cell suspension into flow tube and detect the cells by Flow cytometry.
7. Analyze the result data using FCS Express 7Plus and GraphPad Prism 5 software.

■ Results

The data showed that the expression level of anti-CD19 scFv on the surface of anti-CD19 CAR-293 cells was 100%



2e5 of Anti-CD19 CAR-293 cells were stained with 100 μ L of 1:50 dilution (2 μ L stock solution in 100 μ L FACS buffer) FITC-Labeled Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Cat. No. [FM3-FY45](#)) and isotype control respectively. FITC signal was used to evaluate the binding activity (QC tested).

The protocol can be provided and feel free to send email to cart@acrobiosystems.com to request.

[Click here](#)

Evaluation of Anti-MSLN CAR Expression with PE-labeled MSLN

■ Reagents

PE-labeled Human Mesothelin / MSLN (296-580) Protein (ACROBiosystems, Cat. No. [MSN-HP223](#));

■ Samples


CAR-RC218 cells (Transfected 293 cells expressing the anti-MSLN scFv).

■ Protocol

1. Culture CAR-RC218 cells in DMEM medium with 10% FBS in the CO₂ incubator (at 37°C, 5% CO₂).
2. Harvest the cells and wash the cells once by wash buffer.
3. Count the cells number and the viability, aliquot up 1e6 live cells into each tube.
4. Add 100 µl of diluted PE-labeled Human Mesothelin (296-580) Protein (Cat. No. MSN-HP223) (prepared in dilution buffer at 1:25 dilution) into each tube, incubating at 4°C for 1 hour.
5. Wash the cells 3 times by wash buffer and resuspend the cells in 200 µl PBS per sample.
6. Transfer the cells into flow tube and detect by Flow cytometry.
7. Analyze result using FACS Celesta software and FCS Express 6 Flow software.

■ Results

The data showed that the expression level of anti-MSLN scFv on the surface of R1013-C6 cells was 100 % .



1e6 of the MSLN-CAR-293 cells were stained with 100 µl of 1:25 dilution (4 µl stock solution in 100 µl dilution buffer) of PE- labeled Human Mesothelin / MSLN (296-580) Protein (Cat. No. [MSN-HP223](#)). PE Streptavidin was used as negative control (QC tested).

The protocol can be provided and feel free to send email to cart@acrobiosystems.com to request.

[Click here](#)

V. Product List of CAR-T Targets

Fluorescent-labeled Proteins

Targets	Cat. No.	Product Description
BCMA	<u>BCA-HF254</u>	FITC-labeled Human BCMA, Fc Tag, DMF Filed
BCMA	<u>BCA-HF2H1</u>	FITC-labeled Human BCMA, His Tag, DMF Filed
CD19	<u>CD9-HF2H2</u>	FITC-labeled Human CD19, His Tag, DMF Filed
CD19	<u>CD9-HF251</u>	FITC-labeled Human CD19, Fc Tag, DMF Filed
CD19	<u>CD9-HP2H3</u>	PE-Labeled Human CD19 (20-291) Protein, His Tag
CD22	<u>SI2-HF2H6</u>	FITC-labeled Human Siglec-2 / CD22, His Tag
CD30	<u>CD0-HF2H3</u>	FITC-labeled Human CD30, His Tag
CD4	<u>CD4-HF255</u>	FITC-labeled Human CD4, Fc Tag
CD4	<u>CD4-HF2H7</u>	FITC-labeled Human CD4, His Tag
EGF R	<u>EGR-HF2H5</u>	FITC-labeled EGF R Protein, His Tag
FOLR1	<u>FO1-HF2H8</u>	FITC-labeled FOLR1, His Tag
GPC3	<u>GP3-HF258</u>	FITC-labeled Human GPC3, Fc Tag
GPC3	<u>GP3-HF2H1</u>	FITC-labeled Human GPC3, His Tag
GUCY2C	<u>GUC-HF2H8</u>	FITC-labeled Human GUCY2C / Guanylyl cyclase C Protein, His Tag
Her2	<u>HE2-HF224</u>	FITC-labeled Human Her2, His Tag
Her3	<u>ER3-HF2H5</u>	FITC-labeled Human ErbB3, His Tag
IL13RA2	<u>IL2-HF2H3</u>	FITC-labeled Human IL-13 R alpha 2 Protein, His Tag
MSLN	<u>MSN-HF253</u>	FITC-labeled Human MSLN, Fc Tag
MSLN	<u>MSN-HF223</u>	FITC-Labeled Human MSLN, His Tag
MSLN	<u>MSN-HP223</u>	PE-labeled Human MSLN, His Tag
MSLN	<u>MSN-HG2P4</u>	GFP Fusion Human Mesothelin / MSLN (296-580) Protein, His Tag

Fluorescent-labeled Proteins

Targets	Cat. No.	Product Description
MSLN	<u>MSN-HR2P4</u>	RFP Fusion Human Mesothelin / MSLN (296-580) Protein, His Tag
MSLN	<u>MSN-HF25x</u>	FITC-labeled Human Mesothelin / MSLN (296-580) Protein, Fc Tag
PSMA	<u>PSA-HF244</u>	FITC-labeled Human PSMA, His Tag
SLAMF7	<u>SL7-HF2H7</u>	FITC-labeled Human SLAMF7, His Tag

Biotinylated Proteins

Targets	Cat. No.	Product Description
BCMA	<u>BC7-H82F0</u>	Biotinylated Human BCMA, Fc & Avi Tag, DMF Filed
BCMA	<u>BCA-H82E4</u>	Biotinylated Human BCMA, His & Avi Tag, DMF Filed
CD147	<u>CD7-H82E0</u>	Biotinylated Human EMMPRIN / CD147 Protein, Avi & His Tag
CD19	<u>CD9-H8259</u>	Biotinylated Human CD19, Fc Tag, DMF Filed
CD19	<u>CD9-H82E9</u>	Biotinylated Human CD19 (20-291) Protein, His,Avitag™ (SPR verified)
CD19	<u>CD9-H82F7</u>	Biotinylated Human CD19 (20-291) Protein, Fc,Avitag™
CD22	<u>SI2-H82F8</u>	Biotinylated Human CD22, Fc & Avi Tag
CD30	<u>CD0-H82E6</u>	Biotinylated Human CD30, Avi & His Tag
CD33	<u>CD3-H82E7</u>	Biotinylated Human CD33, Avi & His Tag
CD38	<u>CD8-H82E7</u>	Biotinylated Human CD38, Avi & His Tag
CD4	<u>CD4-H82F3</u>	Biotinylated Human CD4, Fc & Avitag
CD4	<u>CD4-H82E8</u>	Biotinylated Human CD4, His & Avitag
CD70	<u>TN7-H82F4</u>	Biotinylated Human CD70, Avi & Fc Tag
CEACAM5	<u>CE5-H82E0</u>	Biotinylated Human CEACAM-5, His & Avitag
EGF R	<u>EGR-H82E3</u>	Biotinylated Human EGFR, His & Avi Tag

Biotinylated Proteins

Targets	Cat. No.	Product Description
EpCAM	<u>EPM-H8223</u>	Biotinylated Human EpCAM, His Tag
EpCAM	<u>EPM-H8254</u>	Biotinylated Human EpCAM, Fc Tag
EpCAM	<u>EPM-H82E8</u>	Biotinylated Human EpCAM, Avi & His Tag
EpCAM	<u>EPM-H82F9</u>	Biotinylated Human EpCAM, Fc & Avi Tag
FOLR1	<u>FO1-H82E2</u>	Biotinylated Human FOLR1, His & Avi Tag
FOLR1	<u>FO1-H82F9</u>	Biotinylated Human FOLR1, Fc & Avi Tag
GPC3	<u>GP3-H82E5</u>	Biotinylated Human GPC3, His & Avi Tag
Her2	<u>HE2-H822R</u>	Biotinylated Human Her2, His Tag
Her2	<u>HE2-H82E2</u>	Biotinylated Human Her2, His & Avi Tag
Her3	<u>ER3-H8223</u>	Biotinylated Human ErbB3, His Tag
Her3	<u>ER3-H82E6</u>	Biotinylated Human ErbB3, His & Avitag
HGF R	<u>MET-H82E1</u>	Biotinylated Human HGF R, Avi & His Tag
MSLN	<u>MSN-H8223</u>	Biotinylated Human MSLN, His Tag
MSLN	<u>MSN-H826x</u>	Biotinylated Human MSLN, Fc Tag
MSLN	<u>MSN-H82E9</u>	Biotinylated Human MSLN, His & Avi Tag
MSLN	<u>MSN-H82F6</u>	Biotinylated Human MSLN, Fc & Avi Tag
MUC16	<u>CA5-H82F4</u>	Biotinylated Human CA125 / MUC16 Protein, Fc & Avi Tag
Nectin-4	<u>NE4-H82E7</u>	Biotinylated Human Nectin-4, His & Avi Tag
Protein L	<u>RPL-P814R</u>	Biotinylated Recombinant Protein L, His Tag
ROR1	<u>RO1-H82E6</u>	Biotinylated Human ROR1, His & Avi Tag
ROR1	<u>RO1-H82F4</u>	Biotinylated Human ROR1, Fc & Avi Tag
SLAMF7	<u>SL7-H82E0</u>	Biotinylated Human SLAMF7 / CRACC / CD319 Protein, His Tag, Avi Tag
VEGFR2	<u>KDR-H82E5</u>	Biotinylated Human VEGF R2, Avi & His Tag

Unconjugated Proteins

Targets	Cat. No.	Product Description
BCMA	<u>BC7-H5254</u>	Human BCMA, Fc Tag, DMF Filed
BCMA	<u>BCA-H5259</u>	Human BCMA, Llama IgG2b Fc Tag, low endotoxin
BCMA	<u>BCA-H522y</u>	Human BCMA, His Tag, DMF Filed
BCMA	<u>BCA-H5253</u>	Human BCMA, Mouse IgG2a Fc Tag, low endotoxin
CAIX	<u>CA9-H5226</u>	Human Carbonic Anhydrase IX / CA9, His Tag
CD133	<u>CD3-H55H7</u>	Human CD133, His Tag (insect cell)
CD19	<u>CD9-H52H2</u>	Human CD19, His Tag, DMF Filed
CD19	<u>CD9-H5250</u>	Human CD19, Llama IgG2b Fc Tag, low endotoxin
CD19	<u>CD9-H5251</u>	Human CD19, Fc Tag, low endotoxin (Super affinity), DMF Filed
CD19	<u>CD9-H5258</u>	Human CD19 (20-291) Protein, Mouse IgG2a Fc Tag
CD70	<u>TN7-H526x</u>	Human CD27 Ligand / CD70, Fc Tag
CD70	<u>CDL-H525a</u>	Human CD27 Ligand / CD70, Mouse IgG2a Fc Tag, low endotoxin
CD30	<u>CD0-H5229</u>	Human CD30 / TNFRSF8, His Tag
CD30	<u>CD0-H5250</u>	Human CD30 / TNFRSF8, Fc Tag
CD30	<u>TN8-H5250</u>	Human CD30 / TNFRSF8, Llama IgG2b Fc Tag, low endotoxin
CD38	<u>CD8-H5224</u>	Human CD38, His Tag
CD38	<u>CD8-H5255</u>	Human CD38, Fc Tag
CD38	<u>CD8-H5253</u>	Human CD38, Mouse IgG2a Fc Tag, low endotoxin
CD38	<u>CD8-H5252</u>	Human CD38, Llama IgG2b Fc Tag, low endotoxin
CD4	<u>CD4-H5259</u>	Human CD4 (Lys 26 - Pro 396), Fc Tag
CD4	<u>LE3-H5228</u>	Human CD4 (Lys 26 - Trp 390), His Tag
CD7	<u>CD7-H5258</u>	Human CD7, Llama IgG2b Fc Tag

Unconjugated Proteins

Targets	Cat. No.	Product Description
CEACAM5	<u>CE5-H5226</u>	Human CEACAM-5 / CD66e Protein, His Tag
CLL-1	<u>CLA-H5266</u>	Human CLEC12A / MICL / CLL-1 Protein, Fc Tag
CLL-1	<u>CLA-H5245</u>	Human CLEC12A / MICL / CLL-1 Protein, His Tag
EGF R	<u>EGR-H5252</u>	Human EGF R, Fc Tag
EGF R	<u>EGR-H5285</u>	Human EGF R, Strep Tag
EGF R	<u>EGR-H5222</u>	Human EGF R, His Tag
EGF R	<u>EGR-H522a</u>	Human EGF R, His Tag, low endotoxin
CD147	<u>CD7-H5222</u>	Human EMMPRIN / CD147 Protein, His Tag
CD147	<u>CD7-H5259</u>	Human EMMPRIN / CD147 Protein, Fc Tag
EpCAM	<u>EPM-H5223</u>	Human EpCAM / TROP1, His Tag
EpCAM	<u>EPM-H5254</u>	Human EpCAM / TROP1, Fc Tag
FAP	<u>FAP-H5263</u>	Human FAP Protein, Fc Tag
Her3	<u>ER3-H5259</u>	Human ErbB3 / Her3 Protein, Fc Tag
Her3	<u>ER3-H5288</u>	Human ErbB3 / Her3 Protein, Strep Tag
Her3	<u>ER3-H5223</u>	Human ErbB3, His Tag
Her3	<u>ER3-H5257</u>	Human ErbB3 / Her3 Protein, Mouse IgG2a Fc Tag
FOLR1	<u>FO1-H528b</u>	Human FOLR1, Strep Tag
FOLR1	<u>FO1-H5253</u>	Human FOLR1, Fc Tag
FOLR1	<u>FO1-H52H1</u>	Human FOLR1 Protein, His Tag (HPLC-verified)
GPC3	<u>GP3-H5258</u>	Human GPC3, Fc Tag, low endotoxin
GPC3	<u>GP3-H52H4</u>	Human GPC3, His Tag
GUCY2C	<u>GUC-H5257</u>	Human Human GUCY2C, Fc Tag

Unconjugated Proteins

Targets	Cat. No.	Product Description
GUCY2C	<u>GUC-H52H5</u>	Human GUCY2C, His Tag
Her2	<u>HE2-H5287</u>	Human Her2 / ErbB2, Strep Tag
Her2	<u>HE2-H5225</u>	Human Her2 / ErbB2, His Tag
Her2	<u>HE2-H5253</u>	Human Her2 / ErbB2, Fc Tag
HGF R	<u>MET-H5227</u>	Human HGF R / c-MET, His Tag
HGF R	<u>MET-H5256</u>	Human HGF R / c-MET, Fc Tag
IL13RA2	<u>IL2-H5256</u>	Human IL-13 R alpha 2 Protein, Fc Tag (MALS verified)
CD123	<u>ILA-H52H6</u>	Human IL-3 R alpha / CD123, His Tag
CD123	<u>ILA-H5255</u>	Human IL-3 R alpha / CD123, Llama IgG2b Fc Tag, low endotoxin
MSLN	<u>MSN-H5223</u>	Human MSLN (296-580), His Tag
MSLN	<u>MSN-H5253</u>	Human MSLN (296-580), Fc Tag
MSLN	<u>MSN-H522a</u>	Human MSLN (296-580), His Tag, low endotoxin
MUC1	<u>MU1-H5252</u>	Human Mucin-1 / MUC-1 (33-167), Fc Tag
Nectin-4	<u>NE4-H52H3</u>	Human Nectin-4, His Tag
NKG2D	<u>NKD-H5265</u>	Human NKG2D, Fc Tag
PSCA	<u>PSA-H52H6</u>	Human PSCA Protein, His Tag
PSMA	<u>PSA-H52H3</u>	Human PSMA, His Tag
ROR1	<u>RO1-H5222</u>	Human ROR1 (165-305, Frizzled domain), His Tag
ROR1	<u>RO1-H5221</u>	Human ROR1 (39-151, Ig-like domain), His Tag
ROR1	<u>RO1-H5223</u>	Human ROR1 (308-395, Kringle domain), His Tag
ROR1	<u>RO1-H5250</u>	Human ROR1 (30-403), Fc Tag
ROR1	<u>RO1-H522y</u>	Human ROR1 (30-403), His Tag

Unconjugated Proteins

Targets	Cat. No.	Product Description
CD22	<u>SI2-H5228</u>	Human Siglec-2 / CD22 (176-687), His Tag
CD22	<u>CD2-H52H8</u>	Human Siglec-2 / CD22 (20-687), His Tag
CD22	<u>SI2-H525a</u>	Human Siglec-2 / CD22 (20-687), Llama IgG2b Fc Tag, low endotoxin
CD33	<u>CD3-H5226</u>	Human Siglec-3 / CD33, His Tag
CD33	<u>CD3-H5257</u>	Human Siglec-3 / CD33, Fc Tag
CD33	<u>CD3-H5259</u>	Human Siglec-3 / CD33, Llama IgG2b Fc Tag, low endotoxin
SLAMF7	<u>SL7-H5225</u>	Human SLAMF7 / CRACC / CD319 Protein
SLAMF7	<u>SL7-H5256</u>	Human SLAMF7 / CRACC / CD319 Protein, Fc Tag
CD138	<u>SD1-H5228</u>	Human Syndecan-1 / CD138, His Tag
VEGFR2	<u>KDR-H5280</u>	Human VEGF R2 / KDR Protein, Strep Tag
VEGFR2	<u>KDR-H5227</u>	Human VEGF R2, His Tag (HPLC-verified)

Her2 BAFFR LAG-3
Fc Receptor Siglec-10
Biotinylated Protein
PD-L1 VEGF165 CD3 epsilon
CD20 PD-1 BCMA
CD27 PVRIG
CD47 PSMA
CFGL1 TFPI
Siglec-15 Integrin
CD24 CD3E & CD3D CD20
CD19 FcRn PCSK9
IL-2 R alpha
CAR-T Target Protein
Glypican 3 Integrin MICA PD-1
CD30 & CD3G
FcRn ADA Service
EGF R B7-H3 BCMA
Integrin TIGIT TGF-beta 1
Fc 4-1BB Siglec-15
Biotinylated Protein
CD20 CD200 GTR Nectin-4
VEGF165 CD73 FGLI
CD69 Nectin-4
VEGF165 PCSK9 IgG1 Fc CD40 **BCMA** PD-L1
SIRP alpha **ADA Service** PSMA
Nectin-4 Biotinylated Protein CD3E & CD3D Immune Checkpoint Protein
IL-2
SPR /BLI analytical service